



**Gaumard®**  
Simulators for Health Care Education

Leadership Through Innovation®



## Super OB Susie® S500.300

Childbirth Training Torso

- Anatomically accurate pelvic landmarks
- Includes full-term neonate featuring palpable landmarks and articulating limbs
- Lifelike cephalic, breech, shoulder dystocia, operative, and C-section deliveries
- Leopold's, version, PPH management, episiotomy repair, and more



## Super OB Susie® S500.300 | Childbirth Training Torso

### Hands-on obstetric skills training

The new OB SUSIE is a realistic childbirth torso designed to aid educators in teaching labor and delivery management skills to learners of all levels.

- Pelvic examination
- Urinary catheterization
- Cephalic and breech delivery
- Operative vaginal delivery
- Shoulder dystocia management
- Umbilical cord management
- Placenta delivery
- Rectal medications



### Realistic anatomy

The anatomically accurate pelvis and term neonate allows learners to develop clinical skills transferable to scenario-based exercises.

- Adult-sized lower torso with hip-joint articulation
- Uniform skin with realistic elasticity and soft feel
- Pelvic landmarks include ischial spines and coccyx
- Full-term neonate with articulated limbs (18.5 in, 6 lb)
- Neonate landmarks include fontanelle and sutures

### Customizable and adaptable

OB SUSIE's modular design allows you to add new functionalities to meet specific training needs. Choose from a growing list of optional packages.

- PPH management
- Episiotomy repair
- Cervical dilation and effacement trainer
- Delivery Assist Mechanism

### Routine delivery competencies

- Vaginal examination
- Fetal presentation
- Cardinal movements
- Delivery maneuvers
- Cord management
- Placenta delivery



Palpable pelvic landmarks and dilating cervix



Illustrate descent, rotation, and expulsion



Removable abdomen offers an internal view of the birth

### Cord and placenta management

- Nuchal cord
- Cord prolapse
- True knots
- Clamp and cut
- Retained placental fragments



Umbilical cord clamping and cutting



Placenta delivery and examination

## Super OB Susie® S500.300 | Childbirth Training Torso

### OBSTETRIC EMERGENCIES

OB SUSIE can illustrate realistic shoulder dystocia situations to practice advanced delivery maneuvers, teamwork, and communication.

- McRobert's maneuver
- Suprapubic pressure
- Rubin
- Woods' screw
- Arm sweeps




Articulated joints provide realistic resistance



Breech maneuvers: Pinard's, Mauriceau, Ritgen's, Lovset



Forceps indication, application, and traction



Vacuum cup application, suctioning, and traction

### Breech and operative delivery

The neonate's smooth skin and articulated limbs support the use of real instruments and advanced maneuvers.



Soft abdomen and adjustable uterine tone for massage training



Cervix and uterus allow for tamponade placement and hemorrhage control

### PPH management package

Train postpartum care and emergency management skills with the optional PPH package.

- Fundal and bimanual massage
- Uterine tamponade placement
- Urinary catheterization
- Suppository administration



# Super OB Susie® S500.300 | Childbirth Training Torso

## OB SUSIE® Childbirth Training Torso

S500.300 ●●●

The OB SUSIE Childbirth Torso includes everything to get started in one easy-to-use package.



### Neonate Features

- Full-term size and weight (18.5 in, 6 lb.)
- Soft full-body skin
- Nasal and oral cavities
- Palpable fontanelle and sutures
- Articulated limb joints
- Detachable umbilical cord

### Package Contents

- OB SUSIE Torso
- Pregnant abdomen
- Full-term neonate
- Placenta
- (2) Umbilical Cord
- (2) Umbilical Stumps
- (2) Cervices
- (2) Birth Canals
- Mineral Oil
- Talcum Powder
- Blood and urine fill kits

### Torso Features

- Adult-sized lower torso: diaphragm to upper legs
- Smooth and supple skin
- Removable pregnant abdominal cover
- Modular cervix and birth canal
- Bony landmarks including ischial spines and coccyx
- Lifelike placenta with detachable umbilical cord
- Patent urethra with bladder access
- Integrated urine bladder (400 mL)
- Integrated blood reservoir
- Articulating hip-joints
- Patent rectum

## PPH Management Package

S500.300.117 ●●●

The PPH Management Package enables the training of postpartum care and emergency management skills.



### Features

- Realistic post-delivery perineum
- Uterine tamponade placement
- Adjustable uterine tone
- Urinary bladder catheterization
- Suppository administration

### Package Contents

- PPH Perineum
- Boggy Uterus
- PPH Cervix
- Blood concentrate

## Delivery Assist Mechanism

S500.300.929 ●●●

The Delivery Assist Mechanism manual lever makes controlling descent and rotation virtually effortless.



### Features

- Controlled descent
- Fetal rotation control
- Hold during descent

### Package Contents

- Birthing mechanism
- Rear panel
- Fasteners set

## Postpartum Episiotomy Suture Trainer

S500.300.133 ●●●

The episiotomy repair trainers simulate human tissue and support the use of real sutures and instruments.



### Features

- Realistic soft skin
- Supports real sutures
- Durable self-healing skin

### Package Contents

- Mediolateral left episiotomy perineum
- Midline episiotomy perineum

## Request a Quote

[www.gaumard.com/quote](http://www.gaumard.com/quote)

[sales@gaumard.com](mailto:sales@gaumard.com)

Toll-Free USA & Canada

1.800.882.6655

Worldwide 305.971.3790



**Gaumard®**  
Simulators for Health Care Education

Product design and price subject to change without notice. Patented; other patents pending. All trademarks and/or copyright materials are the property of their respective owners. © 2020 Gaumard Scientific. All Rights Reserved. MADE IN USA. 11190027B