

Features	S200 SUSIE SIMON®	S201 SUSIE SIMON® with Ostomy	S205 Simple SIMON™	S206 Simple SUSIE®
General Patient Care				
Bathing and bandaging activity	•	•	•	•
Full body	•	•	•	•
Interchangeable genitalia	•	•	-	-
Eyes open and close	•	•	•	•
Realistic eyes for ophthalmic exercises	•	•	•	•
Realistic urethral passage and bladder for catheterization exercises	•	•	-	-
Upper and lower dentures for oral hygiene	•	•	•	•
Soft, realistic face skin, hands, feet, fingers, and toes	•	•	•	•
Simulated ear canal for otic drops and irrigation	•	•	•	•
Colostomy and ileostomy stomas to practice irrigation	○	•	-	-
Gastrointestinal procedures and enema administration	•	•	-	-
Stylish wig for haircare exercises and surgical draping	•	•	-	•
Amputation stump	○	○	○	○
Set of two decubitus ulcers	○	○	○	○
Ulcerated foot	○	○	○	○
Articulating head, jaw, elbows, wrists, ankles, and knees	•	•	•	•
Injection Training				
IM injection sites in deltoids, quadriceps, and upper gluteal region	•	•	•	•
Advanced Multipurpose Intravenous Training Arm	○	○	○	○
Advanced Intravenous Training Arm	○	○	○	○
Arterial and Venous Training Arm	○	○	○	○
Heart and Lung Sounds				
Site-Specific Heart and Lung sounds with Virtual Stethoscope	○	○	○	○
Airway				
Tracheotomy placement	•	•	•	•
Nasal and oral tube placement	•	•	-	-
NG and OG tube feeding and gastric suction	•	•	-	-
GYN Training				
Realistic vagina and cervix supports douching and pap smear	•	•	-	•
Other				
Detachable and removable internal tanks	•	•	-	-
Bends and detaches at waist for easy storage	•	•	•	•
Latch provides secure seal between ostomies and internal tanks	-	•	-	-
Carrying bag	○	○	○	○
Instruction manual	•	•	•	•