


US-16

NEW

Intravesical Urine Volume Measurement Phantom

Pocket-sized ultrasound is gaining popularity and starting to be utilized more in daily care for examination. This method is preferred as catheterization may cause urinary infection.




FEATURES

1 | 4 variations of interchangeable inserts feature different patient scenarios.


SKILLS

| Handling and manipulation of transducers on bladder ultrasound
| Scanning the bladder to measure urine volume and findings.


KEY FEATURES



Urine volume 50mL



Urine volume 150mL



Urine volume 300mL



Urinary retention + balloon catheter

DESCRIPTIONS

SPECIFICATIONS

Size: 30 x 26 x 18 cm
12 x 10 x 7 inch

Phantom: 1.5kg / 3.3lbs
Body: 3.2kg / 7.1lbs

MATERIALS

Latex free

RECOMMENDED DEVICES

-

SET INCLUDES

- 1 lower torso manikin
- 4 ultrasound blocks
- 1 storage case
- 1 talcum powder
- 1 instruction manual

REPLACEMENT PARTS

-