

Casualty Care Rescue Randy

Life-like training solutions for real-life situations


A full-body manikin, uniquely designed for enhanced realistic training on the procedures that treat the three most preventable causes of death:

1. Massive bleeding
 2. Tension pneumothorax
 3. Airway obstruction
- Compliant with TCCC/TECC scenarios
 - Rugged, reliable, and hyper-realistic
 - 100% hand-made in the USA

Available Blood Pumping System #BPS2-B
 Allows for all four CCRR wounds to bleed simultaneously, and includes additional wounds



- Casualty Care Rescue Randy Manikin
- Hard transfer case
- Outer skin
- Inner skin
- Foot pump with pressure regulator
- 10 trachea
- Three repairable neck skins
- Silicone repair kit



This manikin is intended to be used as an instructional aid in the teaching and demonstrating of live rescue training. Actual product may vary slightly from photo.



01530 512425



shop@simulaidsttd.co.uk



Simulaidsttd™



Simulaids™

Product Name	Casualty Care Rescue Randy – Powered by Strategic Operations Hyper-Realistic® Technology
Addresses Training Needs/Curricula	TCCC – Tactical Combat Casualty Care for all service members TECC – Tactical Emergency Casualty Care for all EMS/EMT/LEO professionals
Configuration	<ul style="list-style-type: none">▪ Full-body manikin with articulated joints▪ Weight-distributed to mimic the human anatomy▪ Holds three to four litres of blood▪ Pumping system with regulator simulating human blood flow pressure (two to three psi)▪ Manual foot-operated blood pumping system allowing three wounds to bleed simultaneously▪ Optional electric deluxe blood pumping system allowing four wounds to bleed simultaneously▪ Blood powder included (makes five gallons)▪ Three repairable neck skins▪ 10 repairable tracheas▪ One repair kit
Dimensions	Male at 72" tall x 25" wide x 12" deep (182.88 cm x 63.5 cm x 30.48 cm)
Product Ship Weight	122 lbs (two boxes)
Product Code/SKU	149-6075
Consumables	<ul style="list-style-type: none">▪ 149-6079 Trachea, Randy - (sold in packs of 10)▪ 149-6080 Repairable neck skin Randy - (sold in packs of three)▪ 149-6081 Silicone repair kit - (one kit)
Additional Options	<ul style="list-style-type: none">▪ IV and suture sleeve▪ IV veins (pack of 10)



01530 512425



shop@simulaids.co.uk



/simulaidsltd